¿About us?


We are a company with more than 8 years of experience in research and technological development focused on medical applications, we work for the well-being of patients, health personnel and to improve the quality of health care services by working on research, development and innovation (R+D+i) for the development of cutting-edge medical technology of the highest quality on the market.

Our research areas are: Endoscopy, Gynecology, Pathology

¿About us?

In endoscopy we have developed the Innovative technological equipment called Clean-Fast which, thanks to its level of innovation, is patent registered in Colombia and internationally before Ompi (PCT).

The image shows two medical professionals in a clinical setting, wearing scrubs and masks. They are positioned on either side of a large, dark-colored machine with a central digital display. The machine has a circular gauge at the top and a screen showing a timer at '00:00:00'. The scene is dimly lit with a blue tint.


Flexible Endoscope Reprocessing

Clean-Fast
Productivity and safety for patients

Clean-Fast

Is a device for reprocessing endoscopes which incorporates the following functions:


- Leak test
- Clean
- Rinse
- Disinfect and dry endoscope chan

Features

Provide a standardized and validated reprocessing cycle in a closed environment:

- Does not use pump for liquid supply
- Temperature and humidity indicator
- Reminder to change disinfectant
- Low energy consumption
- Includes the drying of the channels
- Reduce reprocessing time by 40%
- Communication interface with touch screen user
- Wi-Fi connection
- Ethernet port


- Document process steps automatically (via management software).
- Reduce reprocessing time.
- Provide robust and reproducible reprocessing.
- Minimize environmental pollution
- Facilitate the work carried out by the staff
- Reduce the risk of damage to endoscopes
- Relate preventative maintenance of endoscopes to reprocessing

Communication interface with touch screen user


1. To start the reprocessing process we must enter with a username and a password which guarantees us to keep a record of the person who carried out the reprocessing.

1. Identification of the endoscope to be reprocessed.
2. Leak test
3. Cleaning, rinsing disinfection and drying
4. Indication of the last date that the disinfectant was changed
5. Date and time of reprocessing
6. person who performed the reprocessing


The Clean-Fast linear system provides a **turbulent flow inside the endoscope** channel to improve the sweep of microorganisms and reduce endoscope irrigation and reprocessing time by 40%

The image shows two individuals in a clinical or laboratory setting, wearing full protective gear including masks, face shields, and hairnets. They are positioned on either side of a central piece of equipment. The person on the right is gesturing towards a screen on the device. The entire scene is bathed in a blue light, creating a sterile and high-tech atmosphere. The central device has a digital display showing a timer at 00:00:00 and some status indicators. The background is a plain, light-colored wall.

After performing the reprocessing of the endoscopes, all the information is stored in the cloud where we can access it from a management software and verify the conditions of

Graphic information of the reprocessing

All reprocessing information is stored in the management software with information such as: series of equipment to be reprocessed, date and time of reprocessing, person who performed the reprocessing, environmental conditions, conditions of the disinfectant

SERIE↕	DISPOSITIVO EN FUGA↕	NOMBRE↕	APELLIDOS↕	FECHA DE REALIZACIÓN↕	TERMINADA↕
EG-530FP ENDOSCOPIO CONVENCIONAL	X	Jose Miguel	Soto Cardona	28/03/2022 14:27	✓
EC-530WL3 COLONOSCOPIO C	X	Jose Miguel	Soto Cardona	28/03/2022 12:48	✓
EC-590ZW/L COLONOSCOPIO M	X	Jose Miguel	Soto Cardona	28/03/2022 12:07	✓
EG-530N NASOENDOSCOPIO	X	Jose Miguel	Soto Cardona	28/03/2022 11:38	✓
EC-590ZW/L COLONOSCOPIO M	X	Jose Miguel	Soto Cardona	28/03/2022 11:12	✓
EG-590ZW ENDOSCOPIO MAGNIFICACION	X	Jose Miguel	Soto Cardona	28/03/2022 10:54	✓
EC-590ZW/L COLONOSCOPIO M	X	Jose Miguel	Soto Cardona	28/03/2022 10:15	✓

The reprocessing information is automatically and unlimitedly stored in the software step by step

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Cerrar

Potential market

As for the endoscope reprocessor market in North America, it is valued at **USD 0.59** billion and is expected to reach **USD 0.88** billion by 2026, with a CAGR of 8.5% between 2021-2026. This market is also segmented by product, end user and country (U.S., Canada and the rest of North America) (Market Data Forecast, 2021).

At the Latin American level, the endoscope reprocessor market is valued at **USD 0.16** billion in 2021, and is expected to reach **USD 0.25 billion** in 2026 with a CAGR of 8.58% between 2021 and 2026. (Market Data Forecast, 2021)

Thanks to the level of innovation of Clean-Fast, we have managed to differentiate our equipment on a large scale, unlike some brands in the market where our great differentiator is that we cover 90% of the steps for the reprocessing of endoscopes that start from the endoscope leaks.

The equipment that exists in the market includes cleaning and disinfection the clean-fast incorporates leak test, cleaning, disinfection, drying and additional management software

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Copyright © 2001 - 2017 Fluke Biomedical